

ICOS has a long experience in the treatment of parts and glassware normally used in the pharmaceutical industry.


This new range of ICOS machines show a series of innovative features that satisfy the strict requirements of the pharmaceutical industry (FDA, cGMP, GAMP 5, etc.).

These machines are identified by the item "ML". The "ML" operating process consists of pre-washing, washing, rinsing, drying and cooling phases.




## FIELD OF APPLICATION

These machines are designed for washing, disinfecting and drying of pharmaceutical products such as glassware, pipettes, machine parts, tanks, glass bottles, pill stamps, carboys, tubes, needles, etc.



WATER RECYCLING LOOP WITH PRESSURE AND CONDUCTIVITY CONTROL



AIR HANDLING GROUP WITH TEMPERATURE CONTROL BY RTD



INTERNAL FINISHING OF WASHING CHAMBERS

## MAIN FEATURES

- Distribution of load in racks/carts and baskets for a better performance results
- Washing and rinsing by means of simultaneous combination of water sprayed through nozzles by recycling pump and direct spray of pressurized WFI water; the suspended particles are skimmed off
- Disinfecting by means of dosing different type detergents
- Drying by HEPA filtered hot air at max. 130 °C
- Chamber and service area can be completely inspected/maintained through loading side (dirty area)
- Chamber doors of hinged type with double locking device
- Chamber doors made in accordance with requirements for clean environments
- Possibility of interfacing unloading door with isolator
- Machine equipped with connections for validation tests
- Various types of racks/carts and baskets for any kind of material
- Possibility of application of automatic loading & unloading of racks/trolleys

## CONSTRUCTION FEATURES

- The machine body and associated safety devices are manufactured and tested in compliance with the standards required by the country of destination. Unless there is a specific request, the machines and systems are manufactured and certified in compliance with the European CE mark
- Execution of machine in compliance with cGMP
- Use of sanitary FDA approved and 3A stamped components (valves, filters, piping, tri-clamp ferrules, pumps, etc.).
- Use of 316L stainless steel with roughness less than 0.5  $\mu\text{m}$  (20 micro inches) on all components in contact with the process (machine body, machine doors, spray nozzles bars, process piping, instruments, valves, racks, internal carts, baskets, etc.)
- Respecting of 3-D dead Legs with the application of Zero-Static (0-Dead legs) components
- Respecting of 3 degree slope on process piping
- All internal corners of the chamber are radiused as per cGMP requirement
- Orbital welding on all piping with tri-clamp ferrules ends
- No. 3 Water supplies of "pharma" type with last rinsing water sprayed directly to the spray nozzles/load
- Support frames, external fascia panels and any other "not-process" contact part made of AISI 304 stainless steel
- Machine body insulated with high density ceramic material and with external facing of stainless steel

### CONTROL SYSTEM

The "ML" series are provided with a control system consisting of a controller (PLC), an operating panel (OP) or programmable computer (PC) and printer. The control system manages the machine and its basic parameters, such as phases time, pump pressure, water pressure, water temperature, air temperature, detergents dosing, water conductivity, water PH, water level control it allows an easy interfacing with the operator and makes it possible to customize the process cycles and the manual operations. In addition, the system produces a complete process documentation and assists the users in the operations of machine validation.

Note: A full description of the control system is contained in the brochure.



PLC: SIEMENS or ALLEN-BRADLEY

### STANDARDS AND NORMS

All ICOS machines are in compliance with the most updated standards in the quality of the materials used, in their accurate finishing, construction and process characteristics, such as:

- FDA Regulation 21 CFR part 210 and 211
- FDA Regulation 21 CFR part 11
- cGMP
- GAMP 5
- CE
- UL
- NEMA



OP or PC: SIEMENS or ALLEN-BRADLEY  
SCADA SYSTEM DESIGNED FOR FDA 21 CFR PART 11



RACKS/CARTS AND BASKETS OF DIFFERENT TYPES WITH POSSIBILITY OF CUSTOMIZED TYPE

## PROCESS DATA

The process consists of combined phases and actions designed to provide maximum process flexibility. The choice of times and phase succession makes it possible to obtain the most appropriate cycle for the characteristics of the material to be processed.

| PHASES                | ACTIONS  |
|-----------------------|--|
| loading of material   | through prearranged rack/cart or basket  |
| pre-washing           | water loading, water heating, water recycling and water spraying through nozzles |
| cleaning              | through different type of detergents   |
| drain                 | water drain plus CIP action and PH measurement (optional)                        |
| washing               | water loading, water heating, water recycling and water spraying through nozzles |
| disinfection          | through different type of detergents   |
| drain                 | water drain plus CIP action and PH measurement (optional)                        |
| 1st rinse             | water loading, water heating, water recycling and water spraying through nozzles |
| drain                 | water drain plus CIP action and PH measurement (optional)                        |
| 2nd rinse             | water loading, water heating, water recycling and water spraying through nozzles |
| conductivity control  | on the water recycling loop  |
| drain                 | water drain plus CIP action and PH measurement (optional)                        |
| drying                | hot HEPA filtered air at max. 130°C  |
| cooling               | fresh HEPA filtered air at approx. 20°C  |
| unloading of material | in clean environment by means of prearranged rack/cart or basket                 |

Note: The above phases refer to a standard program.



## ML SERIES - PARTS AND GLASSWARE WASHERS

| Type   | Capacity<br>(litres / cu. ft.) | Chamber dimensions (mm / inches) |             |            | Overall dimensions (mm / inches) |              |             |
|--------|--------------------------------|----------------------------------|-------------|------------|----------------------------------|--------------|-------------|
|        |                                | Width                            | Height      | Length     | Width                            | Height       | Length      |
| ML 100 | 465 / 16.5                     | 700 / 27.5                       | 700 / 27.5  | 950 / 27.5 | 2150 / 84.5                      | 2600 / 102   | 1164 / 45.5 |
| ML 200 | 630 / 22                       | 700 / 27.5                       | 950 / 27.5  | 950 / 27.5 | 2150 / 84.5                      | 2600 / 102   | 1164 / 45.5 |
| ML 300 | 850 / 30                       | 950 / 27.5                       | 950 / 27.5  | 950 / 27.5 | 2400 / 95.5                      | 2600 / 102   | 1164 / 45.5 |
| ML 400 | 3800 / 134                     | 1300 / 51                        | 1950 / 76.5 | 1500 / 59  | 3000 / 118                       | 3800 / 149.5 | 1790 / 70.5 |

ICOS RESERVES THE RIGHT TO MAKE  
PRODUCT CHANGES WITHOUT PRIOR  
NOTIFICATION



**ICOS Pharma S.p.A.**

Via E. Ellero, 15 - Zona Industriale - 33080 Cusano di Zoppola (PN) - Italy  
Phone +39 0434 5772911 - Fax +39 0434 5772900 - e-mail: sales@icospharma.com