TECNinox®

DITEC
Distillatori a
Multiplo effetto
*Multi-effect
distiller*

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Complete solutions

Thanks to the experience TECNinox has gained in the handling of pure fluids, it is undoubtedly the ideal partner for complete water and steam purification systems, from design to execution, through to storage and distribution to the various user. TECNinox has been active in this market sector for several years now and has built up the experience required to develop and construct water for injection generators able to meet the widest Customer requirements for production potential and for checking and monitoring the final steam quality. Together with the Customer user, TECNinox can assess and develop the best technical solution to guarantee continuity of operation and performance, user-friendliness equipment and cost effective operating costs.

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Principle of operation

Distillation process consists of water steaming, followed by steam condensate, to eliminate the most volatile or non-volatile impurities at water boiling temperature.

Multi-effect distillation process includes several effects in which evaporation and condensate usually occur at the same time.

Feed water evaporates in columns consisting of tube condensers and decontamination chambers (so-called "effects") joined by a flange connection and having different pressure and decreasing temperature.

The first effect is fed by an external heat source; in the subsequent effects this energy is totally reused being energetically self-sufficient.

The number of effects is usually included between three and seven.

An increasing number of effects involves more investment costs reducing operative cost at the same time, thanks to higher percentage of recovered energy compared to energy fed in the first effect.

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
DITEC application

Fine chemical, Biotechnology and Pharmaceutical industries need to produce water for injection.

Our range of multi-effect distillers (DITEC) has been designed for the production of Water for injection (WFI) according to cGMP good manufacturing standards in compliance with the European Pharmacopeia and USP guidelines.

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
1 Coil pre-heating

Conventional systems use separate external exchangers connected via external pipework to the decontamination chamber. Our DITEC models are equipped with an integral coil pre-heater in order to:

- Optimize heat recovery
- Simplify maintenance
- Make the system ergonomically easier for accessibility

2 Falling film feed system

Pre-heated purified feed water is fed to multi-tube heat exchangers producing an instantaneous flash evaporation effect.

3 Water droplet separation

In each column water droplets are efficiently separated according to the combination of shear and differential mass.

4 First-effect steam production

All multi-effect distillers usually produce limited Pure Steam in the first column. Since our models are tailor-made according to specific requirements, we can produce any quantity of Pure Steam while distilling WFI.

This technical solution offers several advantages:

- Reduced investment
- Reduced costs and space for installation
- Reduced validation
- Easy handling
- Pure Steam pressure adjustment

5 Automation and control

Access to the PLC based control system is via a touch screen HMI, intuitive and easy to navigate.

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Engineering

TECNinox water distillers belong to the last generation of multi-effect distillers.

The system stability ensures a reliable and automatic production of pyrogen-free and high-quality sterile "Water for Injection" (USP, EP, JP).

Our computer-aided design with 3D modelling, are constantly updated according to the state of the art heat transfer technology.

Distillation columns comply to ASME International rules (U stamp) and PED rules (for Europe), including other design procedures requested by specific countries.

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Service and support

- From the purified water generator to storage and distribution, TECNinox gives its Customers the advantage of working with a single company that is experienced in the handling of pure fluids
- The drawing up and filling in of validation records allows us to offer the Customer complete packages and comprehensive installation qualification (IQ), operation qualification (OQ) and performance qualification (PQ) programs.
- After installation, TECNinox service engineers will provide training courses for the use of the system and provides technical and maintenance support to meet every specific requirement.

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TECNinox reserves the right to modify its own equipment or technical data without previous notification. The above data have a tolerance of $\pm 5\%$.

**Generation
Storage and
Distribution
System for
Purified Media**

Blood Plasma

**Syrup and
Suspension**

**Cream
Gel
Suppository**

Vaccines

**Cytotoxic
Anticancer**

**Parenteral
Formulation
Compounding**

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Documentation

The multi-effect distiller, like all pressure recipients, will be supplied with the conformity certificate to the 97/23/EC PED/ASME directive and to the harmonized 2006/42/EC machine directive.

The relative pressure, temperature features and fluid groups are reported. An user's manual will be provided with the steam generator reporting all maintenance, pre and post installation safety devices, after-sale service details.

Our personnel will deliver the user's manual during training.

Validation package documents will be provided on request in accordance with the pharmaceutical application standard and consist of:

- FAT protocol
- IQ/OQ protocol