
TECNinox®

GVP

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GVP application

Pharmaceutical, Biotechnology, fine Chemical & Cosmetics Industries, Research Institutes and Universities more and more often demand processes with final sterilisation (S.I.P.) for reactors, preparation vessels, autoclaves, filters, piping, etc.

The pure steam produced is saturated and pyrogen free, in accordance with last American USP directives and European Pharmacopeia.

TECNinox's pure steam generators GVP series are designed and built in accordance with the good manufacturing practice cGMP and with FDA standard approved material suitable for pharmaceutical products.

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Principle of operation

The feed-water, pushed by the feeding pump, circulates through an exchanger (optional) to heat the water, then the water goes into the separator column and to the main exchanger (evaporator) pipes side.

The high temperature of industrial steam (shell side) heat the water to the evaporation temperature, creating a thermo siphon circulation type inside the evaporator and the separator column. The steam develops in the column, its low velocity and the height of the column eliminate any possible impure water droplets entrainment.

The pure steam pressure is maintained at the value set by the algorithm (PID) with direct action on the industrial steam modulating valve, evaporator shell side

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
Utilities

- demineralised feeding water free of silica, chlorine, amine and volatile substances.
- dry saturated industrial steam at 9 bar max
- compressed air, without oil, at 6 bar max
- electric energy 400v / 50 hz, 3F + N + T

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Valvola automatica
Automatic Feed Water Control Valve

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Design and construction features

The steam generator is manufactured in AISI 316L stainless steel, PTFE and elastomer gasket with FDA certificate.

The evaporator and the separator column are insulated with mineral wool and cover with AISI 304 stainless steel sheet.

All the components are assembled in a carbon steel tubular painted with resistant epoxy paint.

All fittings, valves, and components are sanitary type, suitable for a non-stop use of pure steam generator.

The control system is inside an electric cabinet. The industrial PLC and the operator panel placed on the front door, control various operating phases, regulating the pure steam pressure outlet, the water level inside the column, the feeding pump, and the visualisation of current parameters.

MODELLO MODEL	ALTEZZA HEIGHT	LUNGHEZZA LENGHT	PROFONDITÀ DEPHT	PESO WEIGHT
GVP 50	1850	900	600	400
GVP 100	2100	1125	600	410
GVP 200	2370	1255	620	560
GVP 300	2500	1255	620	650
GVP 500	2800	1500	650	720
GVP 750	2900	1500	650	830
GVP 1000	3300	1800	700	1120
GVP 1500	3700	2100	800	1520
GVP 2000	4090	2300	1000	1900
GVP 3000	4500	2620	1125	2770
GVP 4000	4600	2750	1200	3520

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TECNinox reserves the right to modify its own equipment or technical data without previous notification.

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Available options

CONDUCTIVITY

Conductivity of feeding water and condensate pure steam can be constantly monitored and recorded. The monitoring system includes an automatic control of set value with an alarm to stop the generator.

ELECTROPOLISHING

Internal surface of under pressure container can be electropolished with a $Ra=0.3$ micron finishing.

FRAME/CONTROL BOARD

It is available on request to manufacture the control cabinet and the framework in AISI 304 stainless steel.

CONDENSING UNIT

The steam generator can be equipped with condensing unit for production of very high quality distilled water for injection. (WFI)

SAMPLING

The sampling unit is composed by a cooling exchanger in order to condense steam and take a sample.

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Documentation

The steam generator will be delivered with the conformity certificate to the 97/23/CE PED directive and 98/37/CE machine directive. The generator will have an external plate with reported the pressure, temperature features and fluid group.

An user's manual will be provided with the steam generator, it reports all the indications regarding maintenance, safety, pre and post installation, troubleshooting and after sales service.

The user's manual will be delivered by our personnel during the commissioning / training.

The validation package documents will be provided on request, in accordance with the pharmaceutical application standard. It is comprising by:

 **FAT protocol**

 **IQ / OQ protocol**

TECNINOX S.r.l.

Via Emilia, 89/A

43015 Noceto, Loc. Sanguinaro (PR) - ITALY

Tel. +39 0521 825324 - Fax +39 0521 825257

www.techninox.it - info@techninox.it

 **TECNINOX®**